

## PRICE LIST

(PRICES EFFECTIVE MARCH 1, 2008)

- RA-250 - *STANDARD* Model.....\$499  
Accommodates patients 4'10" to 6'2"  
Maximum weight 250 lbs.
- HD-400 - *HEAVY DUTY* Model .....\$599  
Accommodates patients 4'10" to 6'2"  
Maximum weight 400 lbs.
- HDT-300 - *TALL* Model.....\$599  
Accommodates patients 6'1" to 6'10"  
Maximum weight 300 lbs.
- HDT-400 - *TALL* Model.....\$599  
Accommodates patients 6'1" to 6'10"  
Maximum weight 400 lbs.
- RA-P-200 - *PEDIATRIC* Model.....\$499  
Accommodates patients 3'2" to 4'6"  
Maximum weight 200 lbs.
- QA-250 - *QUICK ADJUST* Model.....\$699  
Accommodates patients 5'3" to 6'5"  
Maximum weight 250 lbs.
- ATV-250 - *ALL TERRAIN Vehicle*.....\$639  
Accommodates patients 4'10" to 6'2"  
Maximum weight 250 lbs.
- ATV-400 - *ALL TERRAIN Vehicle*.....\$739  
Accommodates patients 4'10" to 6'2"  
Maximum weight 400 lbs.
- ATV-T-300 - *ALL TERRAIN Vehicle* ....\$739  
Accommodates patients 6'1" to 6'10"  
Maximum weight 300 lbs.
- ATV-T-400 - *ALL TERRAIN Vehicle* ....\$739  
Accommodates patients 6'1" to 6'10"  
Maximum weight 400 lbs.
- SW-500 - *STEERABLE* Model .....\$599  
Accommodates patients 4'9" to 6'4"  
Maximum weight 500 lbs.

## Rental Plan

- Model – RA-250.....\$30 per week  
Model – RA-P-200 .....\$30 per week  
HD and HDT models .....\$35 per week  
ATV models .....\$40 per week  
SW-500 models .....\$35 per week  
– *Plus Shipping and Handling* –



*Prices are subject to change without notice.  
Please confirm price when placing your order.*

**Patent #: 5800317**

**The ROLL-A-BOUT is guaranteed.**

**If you are not completely satisfied,  
just return it within one week for a refund  
of the purchase price.**

**Roll-A-Bout Corporation  
3240 Barratts Chapel Rd.  
Frederica, DE 19946**

*www.roll-a-bout.com*

*E-mail:*

*customerservice@  
roll-a-bout.com*

**TELEPHONE ORDERS:**

**TOLL FREE - 888-736-6151**

**Fax 302-335-5484**

# ROLL-A-BOUT

## Eliminates Crutches



***If You Can't Walk About –  
You Need a ROLL-A-BOUT !***

# STEERABLE ROLL-A-BOUT



**SW-500 Model Roll-A-Bout offers the Stability and Safety of our Original Models but with the added feature of Steerable Front Wheels.**

**WHAT IS THE ROLL-A-BOUT?** The patented design of the **ROLL-A-BOUT** is a modern solution to the “age-old crutch.” We have produced the 21st century’s answer to mobility on ONE FOOT! For those who have lost the ability to walk because of a lower leg injury and have tried the **ROLL-A-BOUT**, the comparison is obvious—the **ROLL-A-BOUT** is the VEHICLE OF CHOICE. Crutches are beastly to manage and are even more uncomfortable to use. They are unstable and therefore dangerous to use.

**IS IT EASY TO USE? ... Oh, Yes!** Users of the **ROLL-A-BOUT** like the ease with which it operates — with its four-wheel stability, its hand-operated brakes, and its easy fold-up design. They no longer have to let their body swing by the armpits ... Users rest their knee on a foam-cushioned knee pad (much more comfortable)! Gone are the sore armpits and the lack of stability due to insufficient surface contact. Patients find the **ROLL-A-BOUT** much easier to use than unstable crutches or a walker. It is as easy to use as a child’s *scooter*! To use the **ROLL-A-BOUT** the patient simply places the knee of their injured leg on the special foam knee pad. A second rear-mounted special foam ankle pad supports the injured lower leg. The patient then places one hand on the handlebar with one finger on the bicycle-type hand brake. To move forward simply release the brake and use the uninjured foot to move forward, like a *scooter*. There is no pain inflicted on the patient since the upper leg is carrying all the weight.

*You can just **ROLL-A-BOUT!!***

**WHAT CONDITIONS CAN IT HELP?** The **ROLL-A-BOUT** is best suited for anyone who has an injury below the knee. This includes Heel Spurs, Diabetic Ulcers, Infections, Fractures, Dislocations, Charcot Foot, Ruptured Tendons, Sprained Ankles, or where a lower leg Surgical Operation was necessary. Patients with Amputations and Arthritis may also benefit from the **ROLL-A-BOUT**.

**Why not use CRUTCHES or a WALKER?** NOW there is a better alternative to crutches and walkers ... the **ROLL-A-BOUT**! Crutches and walkers were recommended in cases where doctors wanted patients to keep weight off the foot or leg. Crutches require patients to maintain balance while the weight is shifted in a swinging motion to move forward. *Sore armpits* are the result of the patient suspending the entire body weight on the crutches. *Sprained* and *sore wrists* are also common complaints from people who use crutches. With the **ROLL-A-BOUT**, the injured foot never touches the ground. With the walker, the action is more like *hopping about*, which places a strain on the good foot. In addition, the patient is required to use extra strength in order to keep the foot from touching the ground, which is exhausting to the patient. A patient who lacks the strength to suspend the leg-weight will put the injured foot on the ground, then the injured foot has to bear the weight and can cause healing problems and additional injury.

**The stability of the Roll-A-Bout eliminates these problems, and thus the patient can have a more Enjoyable and Rapid Recovery ... The way it should be!**

**Why not a WHEELCHAIR?** A wheelchair is designed for patients who do not have the strength to stand up on their own. It is important for patients to realize the sooner they get out of the wheelchair, the sooner they will be on the road to recovery! Body strength deteriorates rapidly when in a wheelchair. An extensive physical therapy program is often required for rehabilitation. Cardiovascular fitness, balance and decreased body strength are also concerns when a patient is in a wheelchair. The **ROLL-A-BOUT** is the perfect device for the patient to use while recuperating from a lower-leg injury. The **ROLL-A-BOUT** will restore the patient’s strength, independence, and mobility, while allowing the patient to perform daily hygiene requirements so they can maintain a good quality of life.

**WHERE CAN I GET A ROLL-A-BOUT?** You can obtain a **ROLL-A-BOUT** by placing an order with the factory. You have the option of RENTING or PURCHASING. Rental units require a deposit. We ship units daily via UPS. Many insurance companies will reimburse the costs of rental/purchase when a physician’s prescription and letter stating that the **ROLL-A-BOUT** is a “medical necessity” is submitted. The deposit is returned immediately upon the return of the unit.